


Diagnostics

Precision Perfected

Our dental instrument is meticulously engineered to deliver unparalleled precision in every procedure. With possible treatment.







Front Surface: Rhodium coated, sharp and distortion free mirror image through reflective metal layer directly on glass surface.

Back Surface: The reflective layer lies on the backside of the mirror glass



Fig. 3 41-101-20 Ø 20 mm

Standard Mouth mirror, Back Surface, flat, 12 pieces



Fig. 4

41-102-22 Ø 22 mm

Fig. 5

41-102-24 Ø 24 mm

Standard Mouth mirror, Back Surface, flat, 12 pieces



Fig. 4

41-103-22 Ø 22 mm

Fig. 5

41-103-24 Ø 24 mm

Mouth mirror, magnifying, Back Surface, 12 pieces



Fig. 4

41-104-22 Ø 22 mm

Fig. 5

41-104-24 Ø 24 mm

RHODIUM Mouth mirror, Front Surface, flat, 12 pieces



Fig. 4

41-105-22 Ø 22 mm

Fig. 5

41-105-24 Ø 24 mm

RHODIUM Mega Mouth mirror, Front Surface, flat, 12 pieces



Fig. 4

41-106-22 Ø 22 mm

Fig. 5

41-106-24 Ø 24 mm

RHODIUM Mega Mouth mirror, Front Surface, flat, double-sided, 6 pieces



41-107-22 Ø 22 mm

Fig. 5

41-107-24 Ø 24 mm

RHODIUM Mouth mirror with lip retractor, Front Surface, flat, double-sided, 12 pieces



Fig. 1

41-108-22 Ø 3 mm

RHODIUM MICRO Mouth mirror, Front Surface, flat, 1 piece



41-110-24

RHODIUM MICRO Mouth mirror, Front Surface, flat, 1 piece

3 x 6 mm





Pack of 12 pieces (Packing specimen)



41-109-24 Ø 5 mm

RHODIUM MICRO Mouth mirror, Front Surface, flat, 1 piece
















US (imperial) thread Stainless Steel

42-108-18 1.8 ml 42-108-22 2.2 ml

US (imperial) thread Brass Metric thread Stainless Steel

42-110-18 1.8 ml 42-110-22 2.2 ml

Metric thread Brass Kemp House, 160 City Road, London, EC1V 2NX, United Kingdom info@tuskuk.com | www.tuskuk.com

Imported and Distributed by Prodento
No 637, 5th main, 17th Cross, Indiranagar 2nd Stage, Bangalore 560 038
Info@prodento.in | +91 87920 02278, 98805 11087 | Customer care no: +91 80 2526 6222

